

WINBACK®

MY BODY IS BACK



THERAPY



SPORT

WINBACK

Mindset!


A winning result!

Winback helps speed up rehabilitation for many injured athletes for a winning return to competition.

Present in **25 countries**, Medical Partner with
over **100 International Sports Teams**
& Associations & Growing



MEDICAL PARTNER OF NATIONAL FEDERATIONS



FEEDBACK FROM CUSTOMERS

In less than 5 years, over 3,000 health and sports professionals have started using WINBACK. The community using WINBACK is best placed to talk about our successes so we have asked them for their thoughts.

A winning attitude - a network - experts - results

*With WINBACK ACADEMY, every day,
hundreds of praticionners contribute to the
development of WINBACK technologies.*

Christophe Buée - WINBACK Founder & CEO



*WINBACK isn't just another device, it's a
genuine manual technique that complements
your everyday practice.*

FREDERIC DELACOUR - PHYSIOTHERAPIST

*Winback therapy is part of my daily life now, it
complements and reinforces my manual work.
Its flexibility means I can adapt it to any type
of treatment, from very delicate rehabilitation
of a hand to the thigh of a high performance
athlete*

AUORE CHAMBAUDU - PHYSIOTHERAPIST





I can treat deep tissue with no pain. My patients heal faster, even with chronic conditions like frozen shoulder or osteoarthritis.

**SERGE MAS - OSTEOPATH
PHYSIOTHERAPIST**




The effective draining using Winback technology appealed to me most. The FixPads allow me to keep my hands free and move as usual but with the additional power of Tecartherapy.

LOIC DECHAIINE - PHYSIOTHERAPIST



3 000
physiotherapists
trained in winback TECAR therapy

15 000
winback treatments
every day in France



WINBACK TECARTHERAPY

WINBACK TECAR THERAPY **DRAMATICALLY IMPROVES** the **RESULTS** of your **massage** by reaching deep-seated conditions that hitherto have been deemed unreachable.

The energy developed by your WINBACK machine combines easy handling and intuitive functionality. You concentrate on your manual techniques and the energy compliments each of your expert movements, providing a treatment that is unparalleled in patient comfort and flexibility.



PHYSIOLOGICAL **BIO** ACCELERATOR

REHABILITATION

2X

FASTER ⁽¹⁾

Provide your patients with a faster, more comfortable and pain free recovery.

① REMOVE PAIN FASTER

WINBACK energy blocks pain related to injuries for over 48 hours. Your patient has a sense of well-being and relaxes, making your massage more effective.

② RESTORE FUNCTION FASTER

WINBACK energy releases tension (*contractures, oedema, fibrosis*) by local revascularization. In just 2 minutes you can effortlessly relieve a deep-seated contracture. Blood flows back into the injured tissue once more and the patient regains their range of movement.

③ HEAL INJURIES FASTER

WINBACK energy accelerates the body's natural self-healing system by encouraging sustainable intra and extra cellular exchanges. The benefits of your treatment will carry on after the end of the session.

You provide **the expertise**,
We provide **the technology** !

(1) Results twice as fast on rotator cuff rehabilitation. WINBACK energy comes from "TECAR" technology a term used in scientific studies since 1995. T.E.C.A.R. Therapy. New tendinitis of the rotator cuff: our experience. G. Sanguedolce, C. Venzia, P. Cataldo, G. Mauro Letizia - Chairman of physical medicine and rehabilitation - University of Palermo. Sept 2009. Frequency used 500 kHz

HOW IT WORKS

WINBACK energy is a high frequency current oscillating between **300 kHz and 1 MHz**. It accelerates "natural regeneration" of biological tissues. Its low intensity energy means it is completely non-invasive and 100% natural for the human body.


HEALING
97%
EFFECTIVE ⁽¹⁾

WINBACK energy encourages intra and extra cellular exchanges and therefore faster tissue healing.

REVASCULARIZATION
6X
FASTER ⁽¹⁾

The natural resistance of tissues transforms WINBACK energy into heat (diathermy). This localized increase in heat forces revascularization of restricted areas (fibrosis, contractures) and frees up joint movement.

(1) In 6 sessions, 97% effective on muscle and joint disorders. Evaluation of the effects of multi-frequency capacitive diathermy treatments. S. Pioletti, A. Soldati, F. Speziale, P. Bonifacci, T. Cuzzani, M. Scacchetti, A. Marsotti, S. Alberti, M. Cagnani, R. Marzovillo, A. Garvalli, G. Poletti. Publication Sport & Medicina. Jan/Feb 2009. Frequencies used 500 KHz, 750 KHz and 1 MHz. (2) Clinical Study of Electric Transfer in Hyperthermia pour lumbago treatment K. TAKAHASHI, N. TSUZUKI, K. ZHONG-SHI - Department of orthopaedic surgery - Saltama Medical Center - J. Phys Ther. Sci 11 :45-51, 1999. Frequency used 650 kHz. (3) T.E.C.A.R. Therapy. Postoperative treatment of femoral fracture. A. TERRANOVA, G. VERMIGLIO, S. ARENA, A. CICCIO, S. DI DIO, M. VERMIGLIO - Vol. 44 - Suppl. 1 to No. 3 EUROPA MEDICOPHYSICA - Oct. 2008. Frequency used 500 kHz



WINBACK TECAR THERAPY can act at different tissue depths. Studies led by WINBACK have established 3 working frequencies to reach superficial pathologies, semi deep or deep conditions.

PAIN RELIEF IMMEDIATE AND SUSTAINABLE

WINBACK energy blocks pain-related nerve impulses for over 48 hours. The patient's pain is relieved immediately and for up to 48 hours post treatment.

WINBACK ENERGY IN YOUR HANDS

1 MHz
SUPERFICIAL ACTION

500 KHz
SEMI DEEP ACTION

300 KHz
DEEP ACTION

Contra-indications: transplants, pacemakers and artificial organs, pregnancy, clotting problems, hypersensitivity to temperature changes, burns, infection, cancer, growth cartilage, hypotension, phlebitis, hypersensitivity to pain.

Precautionary principles: osteosynthesis (do not exceed 50% power).

SUPERFICIAL ACTION / DEEP ACTION

Thanks to its **CET and RET** electrodes, you can transfer WINBACK energy into your patient's body in a way that can be local or global, superficial or deep. **Transfer of Energy (TE)** comes from the association between the High Frequency current and electrodes: **CET and RET (CAR)**, making **TECAR**.



SUPERFICIAL ACTION

SOFT TISSUE = CET concentrates energy in the muscles and other fluid tissues. The action is local, superficial or semi deep (DEEPCET)

MUSCLE, FASCIA, APONEVROSIS,
VISCERAL, LYMPHATIC
AND VASCULAR.

DEEP ACTION

HARD TISSUE = RET concentrates energy in all hard tissues such as joints, tendons and bones. The action is poly-articular and deep.

LIGAMENT, JOINT,
TENDON, PERIOSTEUM,
BONES.



**WINBACK FOLLOWS
YOUR MOVEMENTS**



MANUAL MODE (1.0)



HANDS FREE MODE (2.0)



EXPERT MODE (3.0)




EXPRESS MODE (4.0)



INTELLIGENT MODE (5.0)

WINBACK MANUAL MODE 1.0

1 MOBILE ELECTRODE + 1 FIXED ELECTRODE



MOBILE ELECTRODE (CET OR RET)

You move the **MOBILE** electrode to dramatically improve the effects of your tissue massage. With the **CET** mobile electrode, you focus on local soft tissue. With the **RET** mobile electrode, you distribute the energy deep into hard tissues.



FIXED ELECTRODE (RETURN PLATE)

Your second electrode is **FIXED**, it is placed in contact with the patient and closes the electrical circuit with the mobile electrode.




ÉLECTRODE FIXE

ÉLECTRODE MOBILE

The whole range of massage in depth.

The distance between the two electrodes is important and enables treatment of a wide poly-articular areas like the lower limbs. The patient is undeniably comfortable laying on your massage table, leaving you free to work in depth.



An adhesive fixed electrode replaces the traditional return plate. By sticking to the patient's skin, it enables more dynamic treatment (mobilization) or even treating a patient off the treatment bed. Each fixed electrode is designed for single patient use of up to 5-6 treatments complying with strict hygiene standards.

Seated postural stress: Active work on the deep muscles with proprioception are combined with **WINBACK energy** to be faster and more effective.

Sub-quadriceps adhesions: Your working hand emits Winback energy. This means you can quickly relive tissue adhesions.

With **MyBACK*** your hand replaces the mobile electrode. **WINBACK energy** passes into the patient's body through your hand. Completely non-invasive, this technique offers an exceptional experience where you can feel the effects of the energy directly... like a contracture easing within seconds in your hands.



*MyBack bracers are intended to be used for a treatment session of up to 10 minutes maximum.



NEW MyBACK* BRACERS:

New MyBACK* bracers - a must have for manual therapists. The new MyBACK* bracers leave your hands free and are ultra-comfortable to use.

It takes less than a second to put on your MyBACK* bracers and be ready to treat your patient. Moreover, you no longer need cream on your arms. The combination of your skin and the MyBACK* bracers gives you optimal surface contact.

This hands free, time saving devices is a must have for manual therapists.



<-Treating neck pain using new MyBACK bracers (hands free)*

HANDS FREE MODE

WINBACK

2.0

2 FIXED ELECTRODES



FIXED ELECTRODE (RET)

Placing a **FIXED RET** electrode on the patient leaves both your hands completely free.



FIXED ELECTRODE (CET OR RET)

Placing a **FIXED CET or RET** electrode on the patient leaves both your hands completely free.


WINBACK 2.0 is the "Swiss army knife" of WINBACK energy there are so many ways to use it that we can't describe them all. Here are 3 examples of typical treatments currently being used by your colleagues.

DYNAMIC THERAPY

With a wide range of electrodes, all active techniques are compatible with the resistive mode. We can effectively soften the fibrotic tissue of the iliotibial band by stretching, picking, MTP. The patient can also perform all the gym based recovery exercises at the same time.



ONLY FROM
WINBACK



AUTONOMOUS TREATMENT

WINBACK is expanding electrode usage on the patient. With the WINBACK 2.0, it is now possible to apply two fixed electrodes to the patient, one of which replaces the return plate. These electrodes then emit energy alternately. When one emits, the other receives. This new feature enables transmission of homogeneous energy.

Autonomous treatment in 10/15 min. With its low pulse functions and dynamic adaptation of frequencies, Winback offers autonomous application to relieve patients of all tension.

The application of current across fixed electrodes is focused on the areas to be treated for 10-15 minutes. This new principal enables optimisation and development of new protocols.


ASSISTED TREATMENT

In this setting, the electrodes are fixed and the therapist has nothing in his hands for more treatment freedom. The electrodes support the therapist to enhance his manual action.



WINBACK EXPERT MODE 3.0

2 MOBILE ELECTRODES (1 IN EACH HAND)



MOBILE ELECTRODE (CET OR RET)

CET acts on local soft tissue.
RET distributes the energy deep into hard tissues.



MOBILE ELECTRODE (RETURN PLATE)

The second electrode replaces the **WINBACK 1.0** FIXED electrode for better energy management in your hands.

Switch to expert mode with concentrated and dynamic energy. Each electrode alternately emits **WINBACK** energy offering you a focused and very fast treatment (under 10 minutes) on a closely targeted joint or body part. **WINBACK ENERGY** allows you to reach unreachable conditions with ease.

FROZEN SHOULDER:

You work with two active electrodes. **WINBACK** energy is emitted between both your hands. This application is particularly effective in addressing so-called complex joints like the shoulder. This mobile energy means an injury can be tackled from all sides. The whole range of shoulder movements will be restored quickly with unrivalled comfort.



ONLY FROM
WINBACK



WINBACK technology allows therapists to benefit from complete mobility in their movements and to provide faster treatment to the patient.

TECAR therapy allows the therapist to combine high frequency current technology providing cellular stimulation combined with manual therapy.

Thus, the practitioner doesn't need a pre-programmed tool but intuitive technology adapted to their own working methods.

WINBACK provides innovations to accelerate the results from manual therapy approaches. These new technologies are exclusive to **WINBACK** and can reach all injured tissues in record time.



EXPRESS MODE WINBACK 4.0

ONLY FROM
WINBACK

2 ELECTRODES IN ONE

WINBACK 4.0 by Winback is based on MIX high frequency technology, where the two electrodes in contact with the patient are combined into only 1. MIX offers a concentrated and superficial energy that stimulates an immediate reaction in the tissue.

WINBACK 4.0 Both electrodes have been integrated into a single electrode to concentrate the energy on small areas of the body without requiring a lot of power. With 3 electrodes, this mode has emitter/receiver poles in the electrode that are work in pairs through a programmed algorithm.


This mode offers a manual electrode with several dynamic poles that switch alternately between emitting and receiving current. The amount of energy is optimal for cellular stimulation over a surface area of a few centimetres.

6X
FASTER

THAN TRADITIONAL TECAR THERAPY


3X
LESS ENERGY

THAN TRADITIONAL TECAR THERAPY



MIX3 (TRIPOLAR), FOR SUPERFICIAL TREATMENTS (BODY & FOCUS)

MIX3 offers a hand electrode with 3 dynamic poles that switch alternately between emitting and receiving current. Two of the 3 poles are active simultaneously. There are two sizes of **MIX3**: Body and Focus, depending on the size of the treatment area or how delicate the skin is. **MIX3** is dedicated to the superficial layers of tissue. **MIX3** works with a unique frequency of 1 MHz and is particularly suitable for very superficial stimulation.



MIX3 (TRIPOLAR)

WINBACK INTELLIGENT MODE 5.0

ONLY FROM
WINBACK


3 ELECTRODES DISTRIBUTE THE ENERGY

WINBACK Tecartherapy is even further ahead of the game with the introduction of the new **TECARTHERAPY 5.0**: Three electrodes are in contact with the patient, emitting energy through two mobile or fixed electrodes.


This new function means that balanced energy can be transmitted onto two limbs, so, for example, two legs can be treated at the same time. The new set up also means you can work on your massage keeping both of your hands active... Your treatment will be quicker and more effective.

TWO RET ELECTRODES WORK SIMULTANEOUSLY WITH A FIXED RETURN PLATE

2 FIXED ELECTRODES RET & 1 FIXED ELECTRODE NEUTRAL



2 BRACERS RET & 1 FIXED ELECTRODE NEUTRAL



YOUR WINBACK FEATURES

WINBACK energy has features that make your treatment even more efficient. With less energy, you get effective, lasting results with more comfort for your patients. So what are you waiting for?

WINBACK

The right technology
for my patients

TARGETING SUPERFICIAL TO DEEP TISSUE

DEEPCET: QUICKLY RELIEVES MUSCLE TENSION
RET +: MORE COMFORTABLE FOR THE THERAPIST

With BACK 1 & BACK 3

DIATHERMY FOCUS MAXIMUM POWER

LOW PULSE: INFLAMMATION AND LYMPH NODES
SUPER PULSE: MAXIMUM POWER TO RELIEVE TENSION

With BACK 1 & BACK 3





OPTIMISING TREATMENT TIME DETECTING MODE CHANGES

PRESET CET MODE TO A CONTINUOUS 30%

DETECTING WHEN THE ELECTRODES CHANGE BETWEEN CET AND RET FOR ONE MINUTE

TIC FUNCTION: INCREASING OR REDUCING THE INTENSITY OF THE ELECTRODE'S PULSE CONTACT WITH THE PATIENT BY 10%.

100% compatible with the electrodes

DYNAMIC FREQUENCY DYNAMIC AUTOSCAN

DYNAMIC: SCANS ALL TISSUE LAYERS

SWAP: POLY-TISSUE ACTION

BACK 3 only

DOUBLE ACTION TECAR + ELECTROTHERAPY

LOW BEAT: RELIEVE TENSION + CIRCULATION

SUPER BEAT: SUPERFICIAL AREAS

BACK 3 only



THE WINBACK ACADEMY CONNECTS ALL WINBACK USERS TO SHARE & DEVELOP THEIR EXPERTISE

The **WINBACK** academy is committed to sharing its knowledge with you and your colleagues all over the world. Its training courses, forums and social networks are available to support you. You will be able to exchange, learn, and share experiences of **WINBACK** energy with your colleagues who already have it.

3 TRAINING LEVELS TO INCREASE YOUR EXPERTISE:

My START allows you to understand **WINBACK** energy and how to use it with your own techniques. All innovations are presented to help you move forward with your treatments.

My PRACTICE covers all applications of **WINBACK** energy. The 'practical training' places you in a work situation with conditions you treat in your practice. This training allows you to discover new therapeutic domains in rehabilitation, health, sport, and beauty.

My EXPERT develops techniques suited to an area of application. This training is available when you have reached a good level of experience with **WINBACK** energy and considers protocols that integrate other technologies or manual methods.



A COMMUNITY OF 1500 EXPERTS AT YOUR FINGERTIPS

Got a question? Discuss and debate with your colleagues live on our Closed Facebook Community page.

Share your questions, your experiences, your results - get your colleagues' opinions right away.

We've developed the **WINBACK** products to adapt to your movements, your treatment time and your requirements.

THEORETICAL & PRATICAL TRAINING

Your WINBACK training covers the various theoretical and practical aspects of WINBACK Tecartherapy.



TRAINING MANUAL



PROTOCOLS



TRAINING CERTIFICATE



THE BENEFITS OF WINBACK

We've developed the **WINBACK** products to adapt to your movements, your treatment time and your requirements.

INTUITIVE COMPLIMENTS YOUR MANUAL EXPERTISE

WINBACK's energy will reinforce your manual treatment approaches because **it can be used with all known massage techniques:** *manual therapy, structural therapy, osteopathy, vertebral techniques, visceral manipulation, fascia therapy, Chinese medicine, cupping therapy, picking, scar massage, manual lymphatic drainage, posturology, trigger point massage, the Mézières method, muscular chain therapy.*



VERSATILE WITH YOUR NEEDS

Whatever the area (Physio, Sport), there is a tried and tested WINBACK protocol.



THERAPY



SPORT





COMPATIBLE WITH YOUR TREATMENT EQUIPMENT

WINBACK energy **can be used with all of your other treatment equipment** before, after or simultaneously to improve its efficiency. All of your equipment can be used with it: *electrotherapy, cryotherapy, pressure therapy, vibration therapy, mechanical palpate roll, colour therapy, proprioception cushion, Freeman plateau, body weight straps, elastic bands, gym ball, taping.*



FLEXIBLE FITS WELL WITH YOUR TREATMENT PLANNING

Your time is precious and your planning is constantly changing. Organize your day any way you want to, WINBACK allows you to arrange your sessions while remaining efficient. You can position the electrodes for assisted and automated treatment in just 1 minute. In just a few minutes, you can focus the energy locally. Hands-free mode, manual mode, ... It's up to you.



OUR RANGE

THERAPY - SPORT







WINBACK®

PHYSIO

WINBACK energy is compatible with all chronic and acute conditions. In this area, treatment provides visible results as early as the first session (increased range of motion, analgesia, revascularization). The organic return optimizes the healing process while limiting fibrosis.

INDICATIONS

Analgesic and anti-inflammatory, break the vicious circle of painful inflammation / stiffness, free up the range of motion, fibrinolytic action, reduction of oedema, drainage of haematic effusion, improvement of the healing process.

RHEUMATOLOGY

Osteoarthritic conditions: pain relief, deep drainage of oedema, free up range of motion. Chronic tendinopathy, epicondylalgia, neck pain, acute and chronic low back pain, capsulitis, joint stiffness.

PAIN

The rapid elimination of pain is crucial for more effective rehabilitation. The therapist can then start treatment with a faster and deeper action. The analgesic effect lasts for up to 48 hours.

LYMPHATIC DRAINAGE

WINBACK energy optimizes aspiration, pressure and resorption movements and accelerates the manual techniques used in **lymphatic drainage**. Although the effects are immediate, treatment focuses on the mechanical effect, i.e. the rapid and effective removal of excess fluid.

Duration: 10 to 25 minutes - Number: 3-7 sessions
Frequency: 1 per week - Intense: 3 per week

Guy Everts, Belgian physiotherapist: « For five months, I treated a group of 10 patients suffering from fibromyalgia with WINBACK therapy. Today: 4 out of 10 patients I only see every 6 weeks, 4 patients every month and 2 patients every fortnight. My patients generally feel less tired and stronger mentally. »

WINBACK®

S P O R T

WINBACK Sport offers a full range of treatments provided by qualified health professionals in a friendly environment. Sports professionals benefit from exclusive and personalized care.

TRAUMATOLOGY

WINBACK SPORT responds to emergency treatment for acute musculoskeletal injuries. **WINBACK** energy combined with immediate treatment reduces the extent of the trauma and provides a quick, efficient return to physical activity. Trauma causes contusions and is treated immediately to restructure the damaged tissue and stop the inflammatory process. Preoperatively, effusions and oedema are reduced - thus facilitating surgery. Postoperatively, therapeutic complications are reduced.

INDICATIONS

Sprains, muscle tears, fractures, pre-and post-operatively on ligamentoplasties, prostheses and osteosynthesis. Acute post-traumatic (J1-J3 use exclusively non-thermal): analgesic and anti-inflammatory, anti-oedematous. Chronic or sequelae post-traumatic: fibrinolytic action.

RECOVERY

After a competition, a deep tissue massage with **WINBACK** energy effectively relaxes the muscles that have been strained by intense effort. These stimulating and relaxing movements restore the biomechanical structure more quickly while improving muscle tone. Muscle spasm, stiffness, new and old muscle injuries: sprains, tears, strains, muscular contusions.

PREVENTION

The sport physiotherapist's task is to identify weaknesses in the patient's biomechanical functions. He combines strengthening and rehabilitation in an injury prevention programme. He then corrects posture to prevent poor movement triggering a weakness.



100 PARTNER SPORTS TEAMS & FEDERATIONS

FOOTBALL CLUBS - FEDERATION FRANCAISE DE FOOTBALL (FRANCE) - FEDERACAO PORTUGUESA DE FUTEBOL (PORTUGAL) - UNION ROYALE BELGE DES SOCIÉTÉS BELGES (BELGIUM) - OLYMPIC FOOTBALL TEAM (GREECE) - POLAND FOOTBALL FEDERATION (POLAND) - IRAN FOOTBALL FEDERATION - ASOCIACION PARAGUAY FUTBAL (PARAGUAY) **FRANCE** CHAMOIS NORTAIS - ESTAC TROYES - FC BORDEAUX - FC LORIENT - NÎMES OLYMPIQUE - FC RENNES - MHSC MONTPELLIER - PSG PARIS - OGC NICE - OL LYON - FC NANTES - FC TOURS - PARIS FC - RC LENS - FC VALENCIENNES - SCO ANGERS - RED STAR **BELGIUM** KV OSTENDE - SPORTING CHARLEROI **PORTUGAL** RIO AVE FOOTBALL - SPORTING CLUB DE BRAGA **SPAIN** DEPORTIVO LEGANES - VILLA REAL - LEVANTE VALENCIA **UK** GLASGOW RANGERS - NOTTINGHAM FOREST - NEWCASTLE **MIDDLE EAST** AL HILAL RYAD SAUDIA - AL JAZEERA CLUB EMIRATES - AL NASR CLUB EMIRATES - AL WASL FOOTBALL EMIRATES - USMA ALGER **OTHERS** BOCA JUNIOR (ARGENTINA) - BAYERN MUNICH FOOTBALL (GERMANY) - LEGIA WARSAW (POLAND)

BASKETBALL CLUBS - FRENCH FEDERATION BASKET BALL - SPURS (SAN ANTONIO, USA) - CSK (MOSCOW, RUSSIA) FRANCE CHORALE DE ROANNE (FRANCE) - ORLEANS LOIRET BASKET - BBP PERIGUEUX - BASKET LANDES - SHARK ANTIBES - CENTRE FEDERAL BASKET BALL (INSEP, FRANCE)

RUGBY CLUBS - RCT TOULON - MHR MONTPELLIER - BO BIARRITZ - UBB BORDEAUX - OYONNAX RUGBY - STADE TOULOUSAIN TOULOUSE - USAP PERPIGNAN

HANDBALL & VOLLEY BALL - KOREAN VOLLEY BALL ASSOCIATION, SWISS VOLLEY FEDERATION, - PSG HAND (PARIS, FRANCE) - PARIS VOLLEY - KNACK RANSTAD VOLLEY (BELGIUM) - ATLANTIS PALLAVOLO BRESCIA (ITALY)

CYCLING TEAMS - FÉDÉRATION FRANÇAISE DE CYCLISME (FRANCE) - AG2R LA MONDIALE (FRANCE) - LA FRANÇAISE DES JEUX (FRANCE) - FORTUNEO VITAL CONCEPT (FRANCE) - LOTTO SUDAL (BELGIUM) - COFINOGA (FRANCE)

NATIONALE SPORTS FEDERATIONS (medical partner) - FEDERATION FRANÇAISE DE JUDO - FÉDÉRATION FRANÇAISE DE NATATION - FEDERACION ESPANOLA DE TRIATLON - RUSSIAN FENCING FEDERATION - FEDERACION ARGENTINA DE YACHTING - KOREAN ARCHERY FEDERATION

INSTITUTIONS - OLYMPIC COMMITTEE OF SOUTH KOREA - CNEA FONT ROMEU - SPORT MILITARY CLINIC (UAE) - HELLENIC INTERNATIONAL OLYMPIC COMMITTEE





Weighing 4 kilos and with its own carrying case, you can take the Back 1S anywhere.

BACK^{1S}

Winback energy, but lighter

The **BACK 1S** offers all the benefits of **WINBACK** technology with minimum weight for maximum efficiency! At just 4 kilos, therapists can take it anywhere. Its carrying case has been specially developed for this purpose.

Fields of application: THERAPY - SPORT.

FEATURES

Power: CET 400 VA - RET 100 Watts
Frequency: 500 KHz & 300 KHz
Function: DYNAMIC - LOW PULSE
SUPER PULSE - TIC - SCAN

Training: MyPRACTICE

Size: 33x33x15 cm

Weight: 4 kg

Optional carrying case

WINBACK 1.0

CET MOBILE: 3 electrodes of different sizes for focused or radial action

RET MOBILE : 3 electrodes of different sizes for small and large areas

FIXED: Return plate and 10 FLEX adhesive plates

OPTION

WINBACK 2.0 : 1 FIXPAD CET/RET/NEUTRAL

WINBACK 3.0 : MOBILE electrode (return)

WINBACK 5.0 : 2 adhesives plates RET, 2 bracers
RET. Double résistive cable (RET)

*2 year equipment warranty - 1 year accessories warranty.
Manufactured to ISO9001 - ISO13485 - KIWA 1984
standards. Certifications : Medical device*

